



danuButton® Compact XS / Compact Set / Premium Set

The danuButton® is a replacement feeding tube for long-term enteral nutrition and use with completely healed stoma. The low-profile feeding tube is easy to conceal under clothing and especially for active, mobile patients. The extensive accessories of the danuButton® Sets enable you to an easy and comfortable replacement of the gastrostomy probe, additional components are no longer necessary.

The decision for the danuButton® is also the decision for more quality of life through an innovative, high-quality product.



COMFORT

- **Large selection of sizes and lengths** for the individual and perfect fit
- **Extensive range of accessories** in the Sets for easy and convenient replacement of the gastrostomy tube – additional components not necessary
- **High comfort** due to soft, flexible silicone
- **Low-profile design** for more mobility and quality of life
- **Easy cleaning** and maintenance
- **Tilted** balloon valve

SAFETY


- Made of **high-quality medical grade silicone** with high clarity
- All sets are **latex-free, PVC-free, DEHP-free, BPA-free**
- **Radiopaque stripe** along the entire shaft
- **Apple-shaped inner retention balloon with high balloon stability** for long usage, secure fixation and to avoid leakage
- **Flat, rounded tip** enables atraumatic insertion
- **Recessed tip** to avoid trauma
- Secure **anti-reflux** valve

The danumed® Extension Sets included in the danuButton® Compact and Premium Sets are made of kink-resistant, biocompatible polyurethane. They are equipped with ENFit® connectors according to DIN EN ISO 80369-3 and danumed's innovative rotating connector, the ENSwivel®.



The danuButton® Compact Set contains:

- danuButton® Gastrostomy Button
- danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- danumed® Introducer
- danumed® Balloon Syringe 10 ml
- optiLube lubricant
- Gauze Pads (4x)



Fr-Size						
	12	14	16	18	20	24
0.8	DA 1042 1208	DA 1042 1408	DA 1042 1608	DA 1042 1808	DA 1042 2008	DA 1042 2408
1.0	DA 1042 1210	DA 1042 1410	DA 1042 1610	DA 1042 1810	DA 1042 2010	DA 1042 2410
1.2	DA 1042 1212	DA 1042 1412	DA 1042 1612	DA 1042 1812	DA 1042 2012	DA 1042 2412
1.5	DA 1042 1215	DA 1042 1415	DA 1042 1615	DA 1042 1815	DA 1042 2015	DA 1042 2415
1.7	DA 1042 1217	DA 1042 1417	DA 1042 1617	DA 1042 1817	DA 1042 2017	DA 1042 2417
2.0	DA 1042 1220	DA 1042 1420	DA 1042 1620	DA 1042 1820	DA 1042 2020	DA 1042 2420
2.3	DA 1042 1223	DA 1042 1423	DA 1042 1623	DA 1042 1823	DA 1042 2023	DA 1042 2423
2.5	DA 1042 1225	DA 1042 1425	DA 1042 1625	DA 1042 1825	DA 1042 2025	DA 1042 2425
2.7	DA 1042 1227	DA 1042 1427	DA 1042 1627	DA 1042 1827	DA 1042 2027	DA 1042 2427
3.0	DA 1042 1230	DA 1042 1430	DA 1042 1630	DA 1042 1830	DA 1042 2030	DA 1042 2430
3.5	DA 1042 1235	DA 1042 1435	DA 1042 1635	DA 1042 1835	DA 1042 2035	DA 1042 2435
4.0	DA 1042 1240	DA 1042 1440	DA 1042 1640	DA 1042 1840	DA 1042 2040	DA 1042 2440
4.5	DA 1042 1245	DA 1042 1445	DA 1042 1645	DA 1042 1845	DA 1042 2045	DA 1042 2445
5.0	DA 1042 1250	DA 1042 1450	DA 1042 1650	DA 1042 1850	DA 1042 2050	DA 1042 2450
5.5		DA 1042 1455				
6.0		DA 1042 1460		DA 1042 1860		
6.5		DA 1042 1465				

1 Set per box, 5 boxes per shipping box.

The danuButton® Compact XS contains:


- danuButton® Gastrostomy Button
- danumed® Introducer
- danumed® Balloon Syringe 10 ml

Fr-Size						
	12	14	16	18	20	24
0.8	DA 1043 1208	DA 1043 1408	DA 1043 1608	DA 1043 1808	DA 1043 2008	DA 1043 2408
1.0	DA 1043 1210	DA 1043 1410	DA 1043 1610	DA 1043 1810	DA 1043 2010	DA 1043 2410
1.2	DA 1043 1212	DA 1043 1412	DA 1043 1612	DA 1043 1812	DA 1043 2012	DA 1043 2412
1.5	DA 1043 1215	DA 1043 1415	DA 1043 1615	DA 1043 1815	DA 1043 2015	DA 1043 2415
1.7	DA 1043 1217	DA 1043 1417	DA 1043 1617	DA 1043 1817	DA 1043 2017	DA 1043 2417
2.0	DA 1043 1220	DA 1043 1420	DA 1043 1620	DA 1043 1820	DA 1043 2020	DA 1043 2420
2.3	DA 1043 1223	DA 1043 1423	DA 1043 1623	DA 1043 1823	DA 1043 2023	DA 1043 2423
2.5	DA 1043 1225	DA 1043 1425	DA 1043 1625	DA 1043 1825	DA 1043 2025	DA 1043 2425
2.7	DA 1043 1227	DA 1043 1427	DA 1043 1627	DA 1043 1827	DA 1043 2027	DA 1043 2427
3.0	DA 1043 1230	DA 1043 1430	DA 1043 1630	DA 1043 1830	DA 1043 2030	DA 1043 2430
3.5	DA 1043 1235	DA 1043 1435	DA 1043 1635	DA 1043 1835	DA 1043 2035	DA 1043 2435
4.0	DA 1043 1240	DA 1043 1440	DA 1043 1640	DA 1043 1840	DA 1043 2040	DA 1043 2440
4.5	DA 1043 1245	DA 1043 1445	DA 1043 1645	DA 1043 1845	DA 1043 2045	DA 1043 2445
5.0	DA 1043 1250	DA 1043 1450	DA 1043 1650	DA 1043 1850	DA 1043 2050	DA 1043 2450
5.5		DA 1043 1455				
6.0		DA 1043 1460		DA 1043 1860		
6.5		DA 1043 1465				

1 Set per box, 5 boxes per shipping box.

The danuButton® Premium Set contains:

- danuButton® Gastrostomy Button
- danumed® Extension Set for Gastrostomy Button, 30 cm, Bolus-ENFit®, ENSwivel®, straight button connector
- danumed® Extension Set for Gastrostomy Button, 30 cm, Y-ENFit®, ENSwivel®, angled button connector
- danumed® Introducer
- danumed® Balloon Syringe 10 ml
- danumed® Enteral Syringe 60 ml, ENFit®
- optiLube lubricant
- Gauze Pads (4x)



Fr-Size						
	12	14	16	18	20	24
0.8	DA 1041 1208	DA 1041 1408	DA 1041 1608	DA 1041 1808	DA 1041 2008	DA 1041 2408
1.0	DA 1041 1210	DA 1041 1410	DA 1041 1610	DA 1041 1810	DA 1041 2010	DA 1041 2410
1.2	DA 1041 1212	DA 1041 1412	DA 1041 1612	DA 1041 1812	DA 1041 2012	DA 1041 2412
1.5	DA 1041 1215	DA 1041 1415	DA 1041 1615	DA 1041 1815	DA 1041 2015	DA 1041 2415
1.7	DA 1041 1217	DA 1041 1417	DA 1041 1617	DA 1041 1817	DA 1041 2017	DA 1041 2417
2.0	DA 1041 1220	DA 1041 1420	DA 1041 1620	DA 1041 1820	DA 1041 2020	DA 1041 2420
2.3	DA 1041 1223	DA 1041 1423	DA 1041 1623	DA 1041 1823	DA 1041 2023	DA 1041 2423
2.5	DA 1041 1225	DA 1041 1425	DA 1041 1625	DA 1041 1825	DA 1041 2025	DA 1041 2425
2.7	DA 1041 1227	DA 1041 1427	DA 1041 1627	DA 1041 1827	DA 1041 2027	DA 1041 2427
3.0	DA 1041 1230	DA 1041 1430	DA 1041 1630	DA 1041 1830	DA 1041 2030	DA 1041 2430
3.5	DA 1041 1235	DA 1041 1435	DA 1041 1635	DA 1041 1835	DA 1041 2035	DA 1041 2435
4.0	DA 1041 1240	DA 1041 1440	DA 1041 1640	DA 1041 1840	DA 1041 2040	DA 1041 2440
4.5	DA 1041 1245	DA 1041 1445	DA 1041 1645	DA 1041 1845	DA 1041 2045	DA 1041 2445
5.0	DA 1041 1250	DA 1041 1450	DA 1041 1650	DA 1041 1850	DA 1041 2050	DA 1041 2450
5.5		DA 1041 1455				
6.0		DA 1041 1460		DA 1041 1860		
6.5		DA 1041 1465				

1 Set per box, 5 boxes per shipping box.



All danuButton® Sets are:

- Latex-free
- PVC-free
- DEHP-free
- BPA-free