



danumed® Enteral Syringes, Single-Patient Use

danumed® Enteral Syringes, Single-Patient Use, with female ENFit® connector according to EN ISO 80369-3

For the administration of enteral feeds, medicines and liquids, for decompression and flushing

Extensive range of sizes for safe use in neonates, children and adults

Cost-efficient solution with usage for up to 14 days in home use



COMFORT

Comfortable handling for the user

- **Washable and reusable up to 14 days**
- **Transparent syringe cylinder** for viewing the content
- **Precisely legible** and wipe-proof scaling
- **Secure hold** with ribbed grip plate
- **Smooth**, even plunger movement
- **Noticeable safety stop** to prevent unintentional removal of the plunger from the syringe cylinder

danumed® Syringe Cap

One size fits all:
Universal syringe cap for all sizes and versions of danumed® Enteral Syringes.



danumed® Straw

Round-tip enteral straw for safe and clean withdrawal of enteral nutrition, medication and liquids.

With ENFIT® connector.

Available in different sizes and lengths.



SAFETY

Excellent safety for the patient

- **ENFit® connector** to reduce risk of tubing misconnections
- **“Low-Dose” Tip** in 1 ml and 3 ml syringes
- High quality biocompatible material **free of latex, PVC, DEHP and BPA**
- **Validated grading** for high measurement accuracy



danumed® Enteral Syringes, Single Patient Use, are individually sterile packed and available in standard sizes of 1 ml, 3 ml, 5 ml, 10 ml, 20 ml and 60 ml. Further sizes are available on request.

Safe products according to the standards:

DIN EN ISO 80369-3, DIN EN ISO 7886-1, DIN EN ISO 20695

danumed® Enteral Syringes, Single-Patient Use

REF	Product description	Size	Pieces/ sales unit	Latex- free	PVC-free	DEHP-free	Bisphenol- A-free
DA 1022 0110	danumed® Enteral Syringe, single-patient use, with centric tip, low-dose tip and ENFit®	1 ml	100	●	●	●	●
DA 1022 0310	danumed® Enteral Syringe, single-patient use, with centric tip, low-dose tip and ENFit®	3 ml	100	●	●	●	●
DA 1022 0500	danumed® Enteral Syringe, single-patient use, with centric tip and ENFit®	5 ml	100	●	●	●	●
DA 1022 1000	danumed® Enteral Syringe, single-patient use, with centric tip and ENFit®	10 ml	100	●	●	●	●
DA 1022 2000	danumed® Enteral Syringe, single-patient use, with centric tip and ENFit®	20 ml	100	●	●	●	●
DA 1022 6001	danumed® Enteral Syringe, single-patient use, with eccentric tip and ENFit®	60 ml	50	●	●	●	●

All products are reusable and individually sterile packed. Further sizes are available on request.

danumed® Syringe Caps and Enteral Straws

REF	Product description	Size/ length	Pieces/ sales unit	Latex- free	PVC-free	DEHP-free	Bisphenol- A-free
DA 1021 0055	danumed® Syringe Caps for Enteral Syringes, ENFit®	one size	100	●	●	●	●
DA 1021 1210	danumed® Enteral Straw for milk, ENFit®	12 Fr/ 100 mm	30	●	●	●	●
DA 1021 1215	danumed® Enteral Straw for milk, ENFit®	12 Fr/ 150 mm	30	●	●	●	●
DA 1021 1604	danumed® Enteral Straw for milk, ENFit®	16 Fr/ 40 mm	30	●	●	●	●
DA 1021 1615	danumed® Enteral Straw for milk, ENFit®	16 Fr/ 150 mm	30	●	●	●	●
DA 1021 0605	danumed® Enteral Straw for medication, ENFit®	6 Fr/ 50 mm	30	●	●	●	●

All products are disposable and individually sterile packed.

Related products to danumed® Enteral Syringes, Single-Patient Use:



danumed®
Syringe Cap



danumed®
Straw



danumed® Extension Set
for Gastrostomy Button



danuButton®
Gastrostomy Button

danumed Medizintechnik GmbH –
Individual solutions are our standard.

We are an innovative, German medical technology company with the highest quality standards, many years of expertise in enteral nutrition and clear customer focus. We also meet individual customer requests with flexibility and reliability. Personal, competent, collaborative.

danumed Medizintechnik GmbH • Stadthof 18 • 93059 Regensburg • Germany
Tel + 49 (0) 941/ 890 588-80 • Fax + 49 (0) 941/ 890 588-81
info@danumed.com • www.danumed.com



INNOVATORS IN ENTERAL CARE