



Global Trust in German Design & Quality



# danumed® FeedBag

## Empty bag for enteral nutrition and fluid balancing

The danumed® Feedbag is an empty feeding bag with a large filling volume of 1500 ml.

The ENPlus port on the outlet tube allows connection to the danumed® Gravity Sets and to all other commercially available enteral gravity sets with ENPlus connection.



## SAFETY & COMFORT

- **ENPlus port** for connection to enteral feeding systems
- **Large bag volume** of 1500 ml
- **Inscribable bag material** for clear patient identification
- **Scaling** in 100 ml-steps
- With tube clamp for **quick application stop**
- Large topfiller for **convenient filling** of the bag

### danumed® FeedBag

The high-quality, economical version made of latex-free, DEHP- and BPA-free bag material.

REF	Product description	Pieces/sales unit
DA 1051 1000	danumed® FeedBag	30 pieces

The danumed® FeedBag is individually sterile packed. 30 pieces in a shipping box.

The danumed® FeedBags are

- Latex-free
- DEHP-free
- BPA-free

### danumed® FeedBag Premium

The first-class PVC-free version made of EVA bag material and with premium topfiller.

REF	Product description	Pieces/sales unit
DA 1051 1001	danumed® FeedBag Premium	30 pieces

The danumed® FeedBag Premium is individually sterile packed. 30 pieces in a shipping box.

The danumed® FeedBags Premium are

- Latex-free
- PVC-free
- DEHP-free
- BPA-free